

# SQUID as a shot noise thermometer.

September 16, 2015

## Abstract

A Squid at the end of a transmission line is used as a Shot Noise Thermometry. This is done by applying a DC current larger than the critical current across the squid. For larger currents, the resistance of the squid becomes linear. The squid can then be used as a shot noise sample to calibrate the gain and the noise temperature of the circuit.

$$S_I(f, V, T, R) = \frac{k_B T}{R} \left[ \frac{eV + hf}{k_B T} \coth \left( \frac{eV + hf}{2k_B T} \right) + \frac{eV - hf}{k_B T} \coth \left( \frac{eV - hf}{2k_B T} \right) \right] \quad (1)$$

For a digitizer with a bandwidth B and a gain G the Power is equal to:

$$P(f, V, T, R) = GB \left[ k_B T n(f) + \left( \frac{R}{Z_0 + R} \right)^2 S_I(f, V, T, R) \right] \quad (2)$$

$$Tn(f) = \left[ T_{nc}(1 - |T|^2) + \frac{T_{nw}}{G_L} \right] \quad (3)$$

$$T = \frac{Z_0 - R}{Z_0 + R} \quad (4)$$

Where, the squid resistance is R,  $Z_0$  is the impedance of the lines ( $50 \Omega$ ) and Tn is the noise temperature of Network Chain.  $T_{nc}$  is the noise temperature at the first HEMT amplifier and  $T_{nw}$  is the noise temperature of rest of the network.